

# Ki67 indexes in paired cytological smears and cell blocks of canine lymphoma: preliminary results

Filipe Sampaio<sup>1,2</sup>, Carla Marrinhas<sup>1,3</sup>, Luísa Fonte Oliveira<sup>1,4</sup>, Fernanda Malhão<sup>1</sup>, Célia Lopes<sup>1</sup>, Hugo Gregório<sup>4</sup>, Carla Correia-Gomes<sup>5</sup>, Ricardo Marcos<sup>1</sup>, Marta Santos<sup>1</sup>



1 – Cytology and Hematology Diagnostic Services, ICBAS-UP, School of Medicine and Biomedical Sciences, University of Porto, Porto, Portugal; 2 – Laboratório INNO, Braga, Portugal; 3 – Hospital do Baixo Vouga, OneVet Group, Águeda, Portugal; 4 – Anicura CHV Porto, Porto, Portugal; 5 – Animal Health Ireland, Carrick-on-shannon, Ireland



## Introduction


Ki67 is a protein expressed in all cell cycle phases, except G0 and is used as a proliferation marker in canine lymphoma (CL). Ki67 protein can be immunodetected on biopsies and cytology smears of CL. Ki67 index cut-offs with prognostic significance have been defined for histopathology and flow-cytometry. Our goal was to compare Ki67 indexes estimated on previously-stained cytology smears (PSCS) and on matched cell tube blocks (CTB) of CL.

## Methods

- ❖ PSCS and matched CTB of 20 cases CL (18 nodal, 1 intestinal and 1 cutaneous) were retrospectively selected (Fig. 1).
- ❖ Immunocytochemistry was performed using MIB-1 antibody for Ki67 immunolabelling, and Ki67 index was assessed.
- ❖ The Ki67 indexes determined in PSCS and paired CTB were used to classify the CL as low or high proliferative using the cut-off of 12.2% (defined by flow cytometry)<sup>1</sup>.

## Results


- ❖ In 2 PSCS with long time of archive and 1 CTB out 20 no immunolabelling was detected; excluding these cases, Ki67 index median was 30.5% (SD 22.6%) in PSCS and 44.0% (SD 21.7%) in CTB.
- ❖ No significant difference, nor a correlation existed between Ki67 index values on PSCS and matched CTB.
- ❖ Moderate agreement existed between PSCS and CTB for the classification in lymphoma with low vs high proliferation (Ki67 index cut-off 12.2%) – Table 1, Figures 2 and 3.




**Figure 1.** (A) Selection of previously-stained cytology smears and matched cell tube blocks from canine lymphoma cases. (B) Ten fields of each paired case were microphotographed. Ki67 index estimation was performed using a counting tool program by one observer (FS).

**Table 1.** Ki67 index classification (kappa=0.74).

Cytology	Cell-block		Total
	Low	High	
Low	1	4	5
High	2	13	15
<b>Total</b>	<b>3</b>	<b>17</b>	<b>20</b>



**Figure 2.** B-cell lymphomas. Immunolabeling with MIB-1 on previously-stained cytology smears (A and C) and on matched cell tube blocks sections (B and D), respectively. More than 60% of the neoplastic lymphocytes showed positive nuclear reactivity for Ki67 in both cases. Diaminobenzidine chromogen, Haematoxylin counterstain; bar = 20 µm (A and C); and 17 µm (B and D).



**Figure 3.** Two different cases (A and B; C and D) from B-cell lymphomas. Immunolabeling applied to previously-stained cytology smears (A and C) and on matched cell tube blocks sections (B and D), respectively. The first case showed low (7%) and high (45%) positive nuclear reactivity for Ki67 for previously-stained cytology smears and matched cell tube blocks, respectively. The second case showed high (45%) and low (1%) indexes for previously-stained cytology smears and matched cell tube blocks, respectively. These cases represent discrepant positive nuclear reactivities in both specimens, according to the reported cut-off of 12,2%. Diaminobenzidine chromogen, Haematoxylin counterstain; bar = 24 µm (A and C); and 17 µm (B and D).

**Conclusion:** The estimates of the Ki67 index on cytology and on paired CTB in CL cases can vary. However, when a Ki67 index cut-off of 12,2% for defining low vs high proliferative CL is used, the CTB and cytology tend to agree, especially in recent cytology cases.

## Bibliography

1 - Poggi, A., Miniscalco, B., Morello, E., Comazzi, S., Gelain, M.E., Aresu, L. and Riondato, F. (2015), Flow cytometric evaluation of ki67 for the determination of malignancy grade in canine lymphoma. *Vet Comp Oncol*, 13: 475-480. <https://doi.org/10.1111/vco.12078>  
 2 - Poggi A, Miniscalco B, Morello E, Gattino F, Delaude A, Ferrero Poschetto L, Aresu L, Gelain ME, Martini V, Comazzi S, Riondato F. Prognostic significance of Ki67 evaluated by flow cytometry in dogs with high-grade B-cell lymphoma. *Vet Comp Oncol*. 2017 Jun;15(2):431-440. doi: 10.1111/vco.12184. Epub 2016 Jan 21. PMID: 26792094.